

# Disposable



Sold Box/10. Sterile. Single-Use. 10 ft length. White



NAQD1001 \*



NAQD1004 \*



NAQD2201 \*

Monopolar:



NAQD2103 \*



NAQD2407 \*

Bipolar:

- NAQD1001 Monopolar Bovie Plug (male) to Monopolar Pin (female, Straight)
- NAQD1004 Monopolar Bovie Plug (male) to Monopolar Pin (female, 90 Degree)
- NAQD2201 Monopolar Wolf Generator (female) to Wolf Kleppinger (female)
- NAQD2103 Bipolar Valleylab Banana Pin (male) to Storz/Stryker Kleppinger (female)
- NAQD2407 Bipolar Valleylab Banana Pin (male) to Twin Pin (female, 45 Degree)

\*Style available in reusable or disposable

# ELECTROSURGICAL CABLES

Monopolar | Bipolar | Reusable | Disposable

www.capitalmedicalresources.com | p. 614.657.7780 | f. 734.336.4100 | info@capitalmedicalresources.com

# Reusable

Sold Each. Re-usable. Non-Sterile. 10 ft length.



**Monopolar:** NAQI1001 \*

NAQI1003

NAQI1004 \*  
or 52005-76

NAQW2201 \*



**Bipolar:** NAQI2301 \*

NAQI2401

NAQI2402

NAQI2407 \*



NAQI2900

NAQS2101

NAQS2103

NAQW2204

- NAQI1001 Monopolar Bovie Plug (male) to Monopolar Pin (female, Straight)
- NAQI1003 Monopolar Valleylab Banana Pin (male) to Valleylab Banana Pin (male)
- NAQI1004 or 52005-76 Monopolar Bovie Plug (male) to Monopolar Pin (female, 90 Degree)
- NAQW2201 Monopolar Wolf Generator (female) to Wolf Kleppinger (female)
- NAQI2301 Bipolar Valleylab Banana Pin (male) to Wolf Kleppinger (female)
- NAQI2401 Bipolar Valleylab Banana Pin (male) to Twin Pin (female, Straight)
- NAQI2402 Bipolar Fischer Connector (OMS type) to Twin Pin (female Straight)
- NAQI2407 Bipolar Valleylab Banana Pin (male) to Twin Pin (female, 45 Degree)
- NAQI2900 Bipolar Valleylab Banana Pin (male) to Metz/ Mayo Scissors (female)
- NAQS2101 Bipolar Valleylab Banana Pin (male) to Storz Fischer Connector
- NAQS2103 Bipolar Valleylab Banana Pin (male) to Storz Kleppinger (female)
- NAQW2204 Bipolar Wolf Generator (female) to Twin Pin (female, Straight)
- NAQW2202 Tripolar Valleylab Banana Pin (female) to Wolf Generator (female)



**Tripolar:** NAQW2202

Need something you don't see here? Please ask. We may be able to help!

\* Style available in reusable or disposable